SARS-CoV-2 (NP) Ag HS

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IMMUNOASSAY FOR THE HIGH SENSITIVITY QUANTITATIVE DETERMINATION OF PROTEIN NUCLEOCAPSID (NP) ANTIGENS OF SARS-CoV-2 IN SERUM OR PLASMA

For in vitro diagnostic professional use only

1. INTENDED USE

SARS-CoV-2 (NP) Ag HS is a solid phase enzyme immunoassay for the quantitative determination of SARS-CoV-2 nucleocapsid protein (NP) antigens in human serum or plasma samples. The kit contains sufficient reagents for 96 reactions.

2. SUMMARY

In December 2019, a new coronavirus, now officially called the SARS-CoV-2 coronavirus, was identified in Wuhan in China and was immediately linked to the onset of an associated acute respiratory disease called coronavirus 19 disease (COVID19). COVID-19 signs and symptoms may occur 2 to 14 days after infection, including fever, cough, shortness of breath or difficulty breathing, muscle pain and tiredness. In severe cases, the infection can lead to pneumonia, severe acute respiratory syndrome (SARS), kidney failure and death. The nucleocapsid protein (NP) is the most abundant protein on the helical nucleocapsid of coronaviruses, it envelops the entire genomic RNA. NP also interacts with other viral structural proteins to play important roles when entering host cells and assembling and releasing viral particles. Anti-NP antibodies have been shown to be the first and are the predominant antibodies detectable in patient blood samples after coronavirus infection.

3. PRINCIPLE OF THE TEST

SARS-CoV-2 (NP) Ag HS is based on the ELISA principle and is a two-phase immunoassay kit for incubation. The walls of the microwells are coated with human anti-NP monoclonal antibodies, these antibodies can recognize the NP potentially present in human serum or plasma samples. The sample is dispensed into the wells and after an incubation, the NP antigens potentially present in the sample are captured by the anti-protein NP antibodies immobilized on the plastic while the unbound components are washed off in the subsequent washing phase. Next, a detection solution containing HRPO-conjugated anti-SARS-CoV-2NP monoclonal immunoglobulins. The test is subjected to another incubation. A subsequent washing phase removes everything that is not linked to the immunocomplex present on the well wall. A substrate solution containing 3,3',5,5'-tetramethylbenzidine (TMB) is then added, resulting in the formation of a blue color. The color reaction is stopped by a 2M solution of H₂SO₄, transforming the blue color into yellow, the intensity of this coloring is quantified by an absorbance microplate reader at a wavelength of 450 nm. The intensity of the color is proportional to the amount of anti-NP IgG antibodies captured within the wells.

4. WARNINGS AND PRECAUTIONS

- In compliance with article 1 paragraph 2b European directive 98/79 / EC the use of in vitro medical diagnostic is intended by the manufacturer to secure suitability, performance and safety of the product. Therefore the test procedure, the information, the precautions and warnings in the instructions for use must be strictly followed. The

ELVI021AG SARS-nCoV2 (NP) Ag HS

Info: sales@immunospark.com

use of the kits with analyzers and similar equipment has to be validated. Any change in the purpose, in the project, in the composition or structure and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents that may be caused by these reasons. The manufacturer is not liable for any results by visual analysis of patient samples.

- For in vitro diagnostic use only.
- All the human components have been found non reactive for Anti-HIV-Ab, Anti-HCV-Ab e HBsAg. Nevertheless, all materials should still be considered potentially infectious.
- Do not interchange reagents or strips of different production lots.
- Do not use reagents of other manufacturers together with the reagents of this kit.
- Do not use after the expiration date.
- Use only clean equipment.
- Do not interchange bottle caps.
- Close the bottles immediately after use to avoid evaporation and contamination.
- After first opening and subsequent storage check reagents for microbial contamination before use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense with much precision in the wells.

SARS-CoV-2 (NP) Ag HS can only be used by qualified personnel who are thoroughly familiar with the technical work.

5. COMPONENTS OF THE KIT

5.1. Provided materials

Color code	DESCRIPTION	Q.ty	Units	Symbol
-	Microplate coated with monoclonal Anti-SARS-CoV-2 NP: 12 strips breakable in 8 wells, sealed in aluminium foil bag.	1	Pcs	MT PLATE
● Nero	Conjugate antiNP-HRPO (100X) : 0,12 ml of aqueous solution of antibodies anti-SARS-CoV-2 NP conjugated with HRPO.	1	Vial	CONJ
● Rosso	Calibrator : 2 ng of NP recombinant of SARS-CoV-2 in a buffered base. <i>Lyophilized</i> .	1	Vial	CAL
● Arancione	Reaction buffer : 30 ml of buffer to dilute samples and controls.	1	Vial	DIL
● Celeste	Substrate . 12 ml of 3,3',5,5'-Tetramethylbenzidine (TMB); ready to use.	1	Vial	SUBS TMB
○ Bianco	Wash buffer (10X) : 40 ml of concentrated buffer to wash the wells.	1	Vial	WASH BUF
● Grigio	Stop solution : 12 ml of sulphuric acid, 0.2 mol/l, ready to use.	1	Vial	SOLN
	Instructions for users (<i>MELVIO21AG</i>)	1	Pz	-

Positive control is available for sale separately (*REF. ELVIO21AG-CTR*).

5.2. Tools and materials needed but not provided

- Photometer for microplates with 450/620 nm filters
- Incubator at 37 ° C
- Microplate washer
- Micropipettes with disposable tips (10, 100, 200, 1000 µl)
- Vortex-Mixer
- Disposable test tubes
- Support for test tubes
- Deionized or distilled water.
- Timer.

5.3. Storage and stability of the kit

The reagents should be stored between 2-8 ° C. Do not use reagents after the expiration date. The expiry date is printed on the label of each component and on the outer label of the package.

6. SAMPLES COLLECTION AND STORAGE

Serum or plasma samples can be analyzed with EDTA, sodium citrate or heparin. Highly lipemic, icteric or hemolytic specimens are not recommended. Samples with visible microbial contamination should not be used.

If necessary, use vortex on serum or plasma samples at room temperature to ensure homogeneity. Then centrifuge the samples at 10,000 to 15,000 rpm for 5 minutes before assaying to remove the particulate matter. Do not omit this centrifugation step if the samples are cloudy and contain particles.

6.1. Preparation (dilution) of the samples

It is not needed to predilute the sample. But if necessary, this can be diluted in **reaction buffer 1X**.

7. TEST PROCEDURE

7.1. Reagents preparation

- *Bring all reagents to room temperature (20...25°C) before use.*
- The microplate wells are separable. They contain recombinant SARS-CoV-2 protein adhered to the internal walls. The wells, ready for use, must be stored between 2 ... 8 ° C. Place the unused wells in the pouch with the silica desiccant gel. The product is stable until the expiration date if stored between 2-8 ° C.
- The vial of **Conjugate antiNP-HRPO (100X)** contained 0,12 ml of a concentrated solution of antibodies anti-NP conjugated to HRPO, stabilizers, preservatives. Vortex and centrifuge briefly the solution of **antiNP-HRPO (100X)** and dilute the needed quantity in 1:100 rate with **reaction buffer 1X**. 100 µl of **Conjugate antiNP-HRPO 1X** for each well. Prepare only the needed quantity of **Conjugate antiNP-HRPO 1X**. Bring again to 2-8°C the solution of **Conjugate antiNP-HRPO (100X)** immediately after having used the needed quantity.
- **Calibrator**: reconstitute the lyophilized calibrator with 1 ml of **reaction buffer 1X** reconstitute the to obtain the last point of the curve with concentration with 2000 pg/ml. Leave the reconstituted to stand for 10 minutes and shake gently before making the serial dilutions. Using **reaction buffer 1X** perform the following serial dilutions:

Calibrator volume	Reaction buffer 1X volume	Concentration	Calibration point
2000 pg/ml	-	2000,00 pg/ml	7
250 µl of 2000,00 pg/ml	250 µl	1000,00 pg/ml	6
250 µl of 1000,00 pg/ml	250 µl	500,00 pg/ml	5
250 µl of 500,00 pg/ml	250 µl	250,00 pg/ml	4
250 µl of 250,00 pg/ml	250 µl	125,00 pg/ml	3
250 µl of 125,00 pg/ml	250 µl	62,50 pg/ml	2
250 µl of 62,50 pg/ml	250 µl	31,25 pg/ml	1

Reaction buffer 1X is used as Standard 0 (0 ng/ml).

Reconstituted standards can be stored at -20 ° C for up to one month. Any freeze / thaw cycles are to be avoided.

- The vial of **wash buffer (10X)** contains 40 ml of a concentrated buffer, detergents and preservatives. The content is diluted with deionized or distilled water. Prepare a solution 1X of wash buffer mixing 40 ml of **wash buffer (10X)** with 360 ml of distilled water or deionized water. If precipitates are observed in the bottle of **wash buffer (10X)**, heat the bottle in water to 37 ° C until the precipitates disappear. The diluted wash buffer can be stored at 2-8 ° C for up to one month.
- The vial of **substrate** contains 12 ml of 3,3',5,5'-Tetramethylbenzidine (TMB) and ready-to-use hydrogen peroxide. Store in the dark. The solution is colorless or light blue. If it turns blue, it means it is contaminated and can no longer be used. Once opened, the product is stable until the expiration date if stored between 2-8 ° C.
- The vial of **STOP solution** contains 12 ml of sulphuric acid, 2M (R36/38, S26), ready to use. Once opened, the product is stable until the expiration date if stored between 2-8 ° C.

7.2. Procedure note

It is recommended to pipette all samples and controls within a maximum time of 3 minutes.

7.3. Test procedure

Read the instructions carefully before starting the dosage. To obtain valid results it is essential to follow the instructions exactly. The following procedure has been validated for manual execution. For an execution on automatic instrumentation it is recommended to increase the number of washes by 3 to 5 times and to properly evaluate the volume of the solution to be used to wash a sufficiently large surface of the walls of the well and thus avoid interference. First establish the distribution and identification plan of the samples and controls on a worksheet. Insert the required wells into the microplate holder.

Perform the test in the order established by the instructions, without pauses.

Use new, clean tips for each sample and check for manual application. If an automatic processor is used, refer to the instrumentation manual for measures against contamination and cross-reaction between samples.

1. Pipette **100 µl** of calibrators and samples to be tested.
2. Cover the wells with adhesive film.
3. **Incubate 1 hour at room temperature.**
4. At the end of the incubation, remove the film and aspirate the liquid from the wells. Then wash the wells 4 times with **300 µl** of **Wash buffer**. Prevent the solution from overflowing from the wells. The interval between washing and suction must be at least 5 sec. After washing, gently tap the wells with the opening downwards on an absorbent paper to completely remove the liquid.

Warning: Washing is a critical phase. An inaccurate washing determines one

poor test accuracy and a distorted increase in optical densities.

5. Pipette **100µl** of **Conjugate antiNP-HRPO 1X** in each well. Cover the wells with adhesive film.
6. **Incubate 1 hour at room temperature.**
7. Repeat the step 4.
8. Pipette **100µl** of **Substrate** in each well.
9. **Incubate exactly for 15 min at room temperature (20°...25°C) in the dark.**
10. Pipette **100µl** of **Stop solution** in each well, in the same order as the TMB solution. *During incubation the color changes from blue to yellow.*

Warning: *Samples with a very high positive result can cause precipitates dark chromogen! These precipitates influence the reading of the optical densities. IS recommended to dilute the samples one more time as described in 7.1 and to repeat the test.*

For measurement, adjust the photometer for the microplates **to zero** using the substrate-blankco (blank) **in A1**. *If, for technical reasons, it is not possible to adjust the photometer, subtract the absorbance of the white substrate from all the values of the other absorbances.*

Measure the absorbance of all the wells at **450 nm** and enter all the measured values in the worksheet. Where permitted by the functionality of the reading instrument, also use the contribution of the absorbance read at **405 nm** to increase the linearity of the measurement. It is recommended to make a measurement of the dual wavelength optical densities using **620 nm** as the reference length.

Where they were measured in replicate, calculate the average of the absorbances.

8. QUALITY CONTROL

We recommend that you use control samples according to state and federal regulations. The use of control samples is recommended to ensure the validity of the results day by day.

The test must be performed exactly according to the manufacturer's instructions for use. In addition, the user must strictly adhere to the rules of the GLP (Good Laboratory Practice) or other federal, state and local rules and / or laws. This is particularly relevant for the use of control reagents. It is important to always include enough checks in the test procedure to validate the accuracy and precision of the test.

Each microplate must be considered separately in the calculation and interpretation of the assay results, regardless of the number of plates processed simultaneously. The results are calculated by relating the absorbance value of each sample to the Cut-off value.

9. RESULTS CALCULATION

9.1. Test validation

The test can be validated using a positive control supplied separately from Immunospark with REF *ELVI021AG-CTR*.

9.2. Results calculation

Calculate the average of the results of the replicates (if carried out) for calibration points, possible checks and samples.


Build a calibration curve by reporting the average absorbances for each calibration point on the Y axis at the concentrations of the respective standards. It is suggested to trace the curve using a 4-parameter logistic algorithm (4PL) to obtain the best fitting.

9.3 Attended values

The following table shows the optical density values typically obtained from samples tested with ELVI021AG:

Samples	OD450	Concrtations measured of NP (pg/ml)	Concrtations actual of NP (pg/ml)
Siero indiluito di soggetto sano	0.112	-	-
Siero diluito 1:10 di soggetto sano	0.126	-	-
Siero diluito 1:10 con 2000 pg/ml di NP	2.974	1831	1831
Siero diluito 1:4 con COVID19 severa n°1	0.732	418.267	1673.067
Siero diluito 1:4 con COVID19 severa n°2	0.478	248.933	995.733
Siero diluito 1:4 con COVID19 severa n°3	0.432	218.267	873.067
Siero diluito 1:4 con COVID19 severa n°4	0.694	392.933	1571.733
Siero diluito 1:4 con COVID19 severa n°5	0.513	272.267	1089.067

Below is an example of a typical curve resulting from the test:

Warning: this representation should not be used as an evaluation criterion for measurements in practice!!

10. LIMITATIONS OF THE TEST

Contamination by microorganisms or repeated freeze-thaw cycles can alter the absorbance values. Diagnosis of an infectious disease should not be based on a single test result. It is also important to consider the patient's medical history and symptoms. Test results from immunosuppressed patients and infants have limited value.

11. PERFORMANCE OF THE TEST

11.2 Analytical sensitivity

The limit of detection (LOD) of Nucleoprotein of SARS-CoV-2 was determined by adding 3 standard deviations to the average of the optical densities of 20 replicates of the ZERO standard and calculating the corresponding concentrations.

The found value of LOD has been lower than **5 pg/ml**.

11.2 Diagnostic sensitivity

Diagnostic sensitivity is the probability of the test to provide a positive result in the presence of antigens. Diagnostic sensitivity was found **>95%**.

11.2. Diagnostic specificity

Diagnostic specificity is the probability of the test to provide a negative result in the absence of antigens. The diagnostic specificity is **>93%**.

11.3. Accuracy

To calculate accuracy **Intra-assay** 3 positive sera were tested in the same section by replicating each 20 times. The test reproduced the results in a satisfactory way, attesting to the CV of the **5%**.

To calculate accuracy **Inter-assay** 3 positive sera were tested in independent sessions replicating 10 times each. The test repeated the results in satisfactory maneira, settling on the CV of the **6%**.

11.4. Cross-reactivity

7 nCoV19 negative samples from patients with a positive diagnosis of Influenza A / B or RSV were tested with ELVI021AG. No positive results for nCoV19 were found.

11.5 Spiking e recovery

Different concentrations of SARS-CoV2 NP were tapped into the samples described below. Average recoveries have turned out:

Sample type	Average % recovery	Range %
Serum	91,50	80,00-103,00
Plasma	92,00	83,00-106,00
Cel culture medium	85,40	80,00-90,80


11.6 Calibration

The assay is calibrated against a highly purified recombinant SARS-CoV2 Nucleoprotein produced under conditions of extreme quality.

12. BIBLIOGRAFIA/REFERENCES

1. World Health Organization (WHO), WHO Statement Regarding Cluster of Pneumonia Cases in Wuhan, China, Beijing: WHO; 9 Jan 2020.
2. World Health Organization (WHO), Report of the WHO-China Joint Mission on Coronavirus Disease 2019 (COVID-19), WHO; 28 Feb 2020.
3. Fan Y, Zhao K, Shi ZL, Zhou P (March 2019). "Bat Coronaviruses in China". *Viruses*. 11 (3): 210. doi:10.3390/v11030210. PMC 6466186. PMID 30832341.
4. De Groot RJ, Baker SC, Baric R, Enjuanes L, Gorbalenya AE, Holmes KV, Perlman S, Poon L, Rottier PJ, Talbot PJ, Woo PC, Ziebuhr J (2011). "Family *Coronaviridae*". In King AM, Lefkowitz E, Adams MJ, Carstens EB, International Committee on Taxonomy of Viruses, International Union of Microbiological Societies. Virology Division (eds.). *Ninth Report of the International Committee on Taxonomy of Viruses*. Oxford: Elsevier. pp. 806–28. ISBN 978-0-12-384684-6.
5. Kahn, Jeffrey; McIntosh, Kenneth (November 2005), "History and recent advances in Coronavirus discovery", *Pediatric Infectious Disease Journal*, 24 (11): s223–s227, doi:10.1097/01.inf.0000188166.17324.60, archived from the original on 2020-02-05.
6. World Health Organization (28 January 2020). "Clinical management of severe acute respiratory infection when novel coronavirus (2019-nCoV) infection is suspected".

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